

EPMA (JXA-8500F)

전자탐침미소분석기

The FEG-EPMA (Field Emission Gun Electron Probe Micro Analyzer) utilizes X-ray spectrometry and allows for high speed, high accuracy qualitative and quantitative in-depth surface analysis as well as area analysis. The FEG-EPMA employs a patented "in-lens" Schottky type field emission electron gun.

WDS (wavelength dispersive X-ray spectrometer), high probe current, and small probe diameter, the FEG-EPMA is capable of extreme elemental analysis of sub-micron areas.


Model

JEOL (JXA-8500F)

Specifications


- Source : Schottky FEG
- Accelerating voltage : 1 ~ 30 kV
- Image resolution : 3 nm (WD : 11 mm, 30 kV)
- Magnification : 40 ~ 300,000X (WD : 11 mm)
- WDS (3ch) and EDS (1ch)
- Type of crystals : TAP, LDE2, LIF, LIFH, PETJ, PETH

Comparison of X-ray peaks by WDS and EDS



In a WDS X-ray peaks do not overlap, but in an EDS X-ray peaks of neighboring elements extensively overlap each other

Comparison of line analysis by WDS and EDS



Carbon mapping image

Line profile by WDS and EDS

Observation of small particles in a corrosion-resisting steel (NbC) to protect from high




Location L5145 Tel.02-958-5540

Applications


- Topography and composition images
- Qualitative & quantitative analysis
- Line analysis
- Map analysis & Phase analysis

Map analysis & phase analysis



Analysis of inter face layers of Cu electrode (fine eutectic layers in a solder)

Qualitative & quantitative analysis



Qualitative analysis of Tungsten silicide (WSi₂)