

# XRF (X-ray fluorescence)

X-선 형광분석기

## Model

ZSX Primus II

## Specifications


- Upper surface analysis type
- X-ray tube type
  - Rh target [end-window type]
- Max. load : 4 kW (60 kV, 150 mA)
- Analyzing crystal
  - LiF200, PET, GE, RX25

## Applications


- Qualitative analysis of unknown sample
- Quantitative analysis
- Semi-Quantitative analysis of unknown sample
- Micro mapping analysis
- Thin film analysis




Location L1134 Tel.02-958-5958/5427



## Quantitative analysis of Mg metal



## Thin film analysis



- Distribution of material in micro focus
- Qualitative & Quantitative analysis

## Micro mapping analysis of Mg alloy

