

## ICP-OES (Inductively Coupled Plasma Optical Emission Spectrometer)

유도 결합 플라즈마 분광분석기

### Model

iCAP 6000 Series (Thermo, US)

### Specifications

- Wavelength range : 166-847 nm
- Fuel : argon gas
- Plasma : axially viewed mode

### Applications

- Environmental : water, air, soil, waste
- Material : ceramic, metal alloy, plastic
- Biological and clinical
- Agriculture and food
- Geochemical and metallurgical



## ICP-OES (Inductively Coupled Plasma Optical Emission Spectrometer)

유도 결합 플라즈마 분광분석기

### Model

ICP OES 720 (Agilent, US)

### Specifications

- Wavelength range : 165-800 nm
- Large format programmable array detector
- Fuel : Argon gas
- Plasma: Axially viewed mode

### Applications

- Environmental : water, air, soil, waste
- Material : ceramic, metal alloy, plastic
- Biological and clinical
- Agriculture and food
- Geochemical and metallurgical

